

For Research Use Only

ORC3L Polyclonal antibody, PBS Only

Catalog Number:32429-1-PBS



Basic Information

Catalog Number: 32429-1-PBS	GenBank Accession Number: NM_012381.3	Purification Method: Antigen affinity Purification
Concentration: 1 mg/ml	GeneID (NCBI): 23595	
Source: Rabbit	UNIPROT ID: Q9UBD5	
Isotype: IgG	Full Name: origin recognition complex, subunit 3-like (yeast)	
Immunogen Catalog Number: AG37481	Calculated MW: 66 kDa, 82 kDa	
	Observed MW: 66 kDa, 82 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human

Background Information

ORC3L, also known as Origin Recognition Complex Subunit 3-Like, is a human gene that encodes a protein which is a subunit of the Origin Recognition Complex (ORC). The ORC is a highly conserved six-subunit protein complex that is essential for the initiation of DNA replication in eukaryotic cells. Studies in yeast have shown that ORC binds specifically to origins of replication and serves as a platform for the assembly of additional initiation factors such as Cdc6 and Mcm proteins.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

