

For Research Use Only

PCNX1 Polyclonal antibody

Catalog Number: 32431-1-AP



Basic Information

Catalog Number: 32431-1-AP	GenBank Accession Number: BC140784	Purification Method: Antigen affinity Purification
Concentration: 350 ug/ml	GeneID (NCBI): 22990	Recommended Dilutions: WB 1:1000-1:4000
Source: Rabbit	UNIPROT ID: Q96RV3	
Isotype: IgG	Full Name: pecanex homolog (Drosophila)	
Immunogen Catalog Number: AG37954	Calculated MW: 259 kDa	
	Observed MW: 240 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : A-253 cells,
Species Specificity: human	

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 32431-1-AP (PCNX1 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.