

For Research Use Only

PCNX1 Polyclonal antibody, PBS Only

Catalog Number: 32431-1-PBS



Basic Information

Catalog Number: 32431-1-PBS	GenBank Accession Number: BC140784	Purification Method: Antigen affinity Purification
Concentration: 1 mg/ml	GeneID (NCBI): 22990	
Source: Rabbit	UNIPROT ID: Q96RV3	
Isotype: IgG	Full Name: pecanex homolog (Drosophila)	
Immunogen Catalog Number: AG37954	Calculated MW: 259 kDa	
	Observed MW: 240 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 32431-1-AP (PCNX1 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 32431-1-PBS in a different storage buffer formulation.