

For Research Use Only

ZCCHC4 Polyclonal antibody

Catalog Number: 32481-1-AP



Basic Information

Catalog Number:

32481-1-AP

Concentration:

500 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34661

GenBank Accession Number:

BC157834

GeneID (NCBI):

29063

UNIPROT ID:

Q9H5U6

Full Name:

zinc finger, CCHC domain containing

4

Observed MW:

55 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB: 1:1000-1:8000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HeLa cells, Jurkat cells, Raji cells, Ramos cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 32481-1-AP (ZCCHC4 antibody) at dilution of 1:4000 incubated at room temperature for 1.5 hours.