

For Research Use Only

ZCCHC4 Polyclonal antibody, PBS Only

Catalog Number: 32481-1-PBS



Basic Information

Catalog Number:

32481-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34661

GenBank Accession Number:

BC157834

GeneID (NCBI):

29063

UNIPROT ID:

Q9H5U6

Full Name:

zinc finger, CCHC domain containing

4

Observed MW:

55 kDa

Purification Method:

Antigen affinity Purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 32481-1-AP (ZCCHC4 antibody) at dilution of 1:4000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 32481-1-PBS in a different storage buffer formulation.