

For Research Use Only

ZNF280B Polyclonal antibody, PBS Only

Catalog Number: 32519-1-PBS



Basic Information

Catalog Number:

32519-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG38544

GenBank Accession Number:

NM_080764.3

GeneID (NCBI):

140883

UNIPROT ID:

Q86YH2

Full Name:

zinc finger protein 280B

Calculated MW:

62kDa, 543aa

Observed MW:

60 kDa

Purification Method:

Antigen affinity Purification

Applications

Tested Applications:

WB, IP, Indirect ELISA

Species Specificity:

human

Background Information

ZNF280B belongs to the family of zinc finger proteins, a transcription factor that is mainly involved in the regulation of gene expression. In prostate cancer cells, ZNF280B induces p53 nuclear export, leading to subsequent proteasomal degradation. Therefore, ZNF280B has pro-growth and pro-survival functions in prostate cancer cells. In addition to this, ZNF280B has a promotional role in gastric and lung cancer development, and its potential as a tumor therapeutic target can be explored in the future.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 32519-1-AP (ZNF280B antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 32519-1-PBS in a different storage buffer formulation.



IP result of anti-ZNF280B (IP:32519-1-AP, 4ug; Detection:32519-1-AP 1:1500) with K-562 cells lysate 2440 ug. This data was developed using the same antibody clone with 32519-1-PBS in a different storage buffer formulation.