

For Research Use Only

# FAM40A Polyclonal antibody, PBS Only

Catalog Number:32554-1-PBS



## Basic Information

Catalog Number:

32554-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG38142

GenBank Accession Number:

BC121793

GeneID (NCBI):

85369

UNIPROT ID:

Q5VSL9

Full Name:

family with sequence similarity 40, member A

Calculated MW:

837 aa, 96 kDa

Observed MW:

95-110 kDa

Purification Method:

Antigen affinity Purification

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, rat

## Background Information

FAM40A plays a role in the regulation of cell morphology and cytoskeletal organization. It is required in the cortical actin filament dynamics and cell shape.

## Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 32554-1-AP (FAM40A antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 32554-1-PBS in a different storage buffer formulation.