

For Research Use Only

ZSCAN2 Polyclonal antibody, PBS Only

Catalog Number: 32581-1-PBS



Basic Information

Catalog Number:

32581-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG38126

GenBank Accession Number:

BC136342

GeneID (NCBI):

54993

UNIPROT ID:

Q7Z7L9

Full Name:

zinc finger and SCAN domain
containing 2

Calculated MW:

70 kDa

Observed MW:

70 kDa

Purification Method:

Antigen affinity Purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Background Information

Zinc finger and SCAN domain-containing protein 2 (ZSCAN2, also known as ZFP29 and ZNF854) plays an important role in the maintenance of the SSC (spermatogonial stem cells) pool in vitro and is integral for robust regeneration of spermatogenesis in vivo following exposure to clastogens. Loss of ZSCAN2 expression alters key cellular processes in undifferentiated spermatogonia, such as translation, chromatin modification, and ubiquitin conjugation (PMID: 35244684).

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 32581-1-AP (ZSCAN2 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 32581-1-PBS in a different storage buffer formulation.