

For Research Use Only

# TAF7L Polyclonal antibody

Catalog Number: 32590-1-AP



## Basic Information

Catalog Number:

32590-1-AP

Concentration:

400 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG38542

GenBank Accession Number:

NM\_024885.3

GeneID (NCBI):

54457

UNIPROT ID:

Q5H9L4

Full Name:

TAF7-like RNA polymerase II, TATA box binding protein (TBP)-associated factor, 50kDa

Calculated MW:

53kDa, 462aa

Observed MW:

53 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB 1:500-1:3000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : HeLa cells, mouse testis tissue

## Background Information

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 32590-1-AP (TAF7L antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.