

For Research Use Only

# TAF7L Polyclonal antibody, PBS Only

Catalog Number: 32590-1-PBS



## Basic Information

Catalog Number:

32590-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG38542

GenBank Accession Number:

NM\_024885.3

GeneID (NCBI):

54457

UNIPROT ID:

Q5H9L4

Full Name:

TAF7-like RNA polymerase II, TATA box binding protein (TBP)-associated factor, 50kDa

Calculated MW:

53kDa, 462aa

Observed MW:

53 kDa

Purification Method:

Antigen affinity Purification

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse

## Background Information

## Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 32590-1-AP (TAF7L antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 32590-1-PBS in a different storage buffer formulation.