

For Research Use Only

HBA2 Polyclonal antibody

Catalog Number: 32637-1-AP



Basic Information

Catalog Number:

32637-1-AP

Concentration:

500 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG37751

GenBank Accession Number:

BC008572

GeneID (NCBI):

3040

UNIPROT ID:

P69905

Full Name:

hemoglobin, alpha 2

Observed MW:

15 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB: 1:5000-1:50000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : mouse liver tissue, rat liver tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 32637-1-AP (HBA2 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.