

For Research Use Only

EID2 Polyclonal antibody

Catalog Number: 32693-1-AP



Basic Information

Catalog Number:

32693-1-AP

Concentration:

450 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG38433

GenBank Accession Number:

NM_153232.3

GeneID (NCBI):

163126

UNIPROT ID:

Q8N6I1

Full Name:

EP300 interacting inhibitor of differentiation 2

Calculated MW:

25kDa, 236aa

Observed MW:

28 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

IP: 0.5-4.0 ug for 1.0-3.0 mg of total protein lysate

Applications

Tested Applications:

IP, ELISA

Species Specificity:

human

Positive Controls:

IP: HEK-293 cells,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



IP result of anti-EID2 (IP:32693-1-AP, 4ug;
Detection:32693-1-AP 1:500) with HEK-293 cells
lysate 1800 ug.