

For Research Use Only

ZNF449 Polyclonal antibody

Catalog Number: 32701-1-AP



Basic Information

Catalog Number:

32701-1-AP

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG36778

GenBank Accession Number:

BC060768

GeneID (NCBI):

203523

UNIPROT ID:

Q6P9G9

Full Name:

zinc finger protein 449

Observed MW:

70 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

IP: 0.5-4.0 ug for 1.0-3.0 mg of total protein lysate

Applications

Tested Applications:

IP, ELISA

Species Specificity:

human

Positive Controls:

IP : HEK-293 cells,

Background Information

ZNF449 (Zinc Finger Protein 449) is a member of the zinc-finger family of transcription factors. It is encoded by a gene located on the X chromosome (Xq26.3) and contains an N-terminal SCAN domain and seven C2H2-type zinc finger motifs at the C-terminus.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



IP result of anti-ZNF449 (IP:32701-1-AP, 4ug;
Detection:32701-1-AP 1:600) with HEK-293 cells
lysate 1280 ug.