

For Research Use Only

ZNF449 Polyclonal antibody, PBS Only

Catalog Number: 32701-1-PBS



Basic Information

Catalog Number:

32701-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG36778

GenBank Accession Number:

BC060768

GeneID (NCBI):

203523

UNIPROT ID:

Q6P9G9

Full Name:

zinc finger protein 449

Observed MW:

70 kDa

Purification Method:

Antigen affinity Purification

Applications

Tested Applications:

IP, Indirect ELISA

Species Specificity:

human

Background Information

ZNF449 (Zinc Finger Protein 449) is a member of the zinc-finger family of transcription factors. It is encoded by a gene located on the X chromosome (Xq26.3) and contains an N-terminal SCAN domain and seven C2H2-type zinc finger motifs at the C-terminus.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



IP result of anti-ZNF449 (IP:32701-1-AP, 4 μ g; Detection:32701-1-AP 1:600) with HEK-293 cells lysate 1280 μ g. This data was developed using the same antibody clone with 32701-1-PBS in a different storage buffer formulation.