

For Research Use Only

FAM171B Polyclonal antibody

Catalog Number: 32704-1-AP



Basic Information

Catalog Number:

32704-1-AP

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG37311

GenBank Accession Number:

BC060872

GeneID (NCBI):

165215

UNIPROT ID:

Q6P995

Full Name:

family with sequence similarity 171, member B

Observed MW:

55 kDa

Purification Method:

Antigen affinity Purification

Applications

Tested Applications:

ELISA

Species Specificity:

human

Background Information

FAM171B belongs to the FAM171 family. It is an evolutionarily conserved protein with a wide range of expression across various tissues. FAM171B may play a role in stimulating immune infiltration and immune responses in diseases such as pulmonary arterial hypertension.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data