

For Research Use Only

# FAM171B Polyclonal antibody, PBS Only

Catalog Number: 32704-1-PBS



## Basic Information

Catalog Number:

32704-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG37311

GenBank Accession Number:

BC060872

GeneID (NCBI):

165215

UNIPROT ID:

Q6P995

Full Name:

family with sequence similarity 171, member B

Observed MW:

55 kDa

Purification Method:

Antigen affinity Purification

## Applications

Tested Applications:

Indirect ELISA

Species Specificity:

human

## Background Information

FAM171B belongs to the FAM171 family. It is an evolutionarily conserved protein with a wide range of expression across various tissues. FAM171B may play a role in stimulating immune infiltration and immune responses in diseases such as pulmonary arterial hypertension.

## Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

---

## Selected Validation Data