

For Research Use Only

TRMT9B Polyclonal antibody

Catalog Number:32716-1-AP



Basic Information

Catalog Number:

32716-1-AP

Concentration:

400 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG36370

GenBank Accession Number:

BC016633

GeneID (NCBI):

57604

UNIPROT ID:

Q9P272

Full Name:

chromosome 8 open reading frame 79

Observed MW:

45 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB: 1:1000-1:8000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HEK-293 cells, HepG2 cells, L02 cells

Background Information

TRMT9B is one of two metazoan homologs of yeast tRNA methyltransferase 9 (Trm9), which methylates tRNA wobble uridines (PMID: 37642556).

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 32716-1-AP (C8orf79 antibody) at dilution of 1:4000 incubated at room temperature for 1.5 hours.