

For Research Use Only

UBXN7 Polyclonal antibody

Catalog Number:32722-1-AP



Basic Information

Catalog Number: 32722-1-AP	GenBank Accession Number: NM_015562.1	Purification Method: Antigen affinity Purification
Concentration: 400 ug/ml	GeneID (NCBI): 26043	Recommended Dilutions: WB: 1:1000-1:8000
Source: Rabbit	UNIPROT ID: O94888	
Isotype: IgG	Full Name: UBX domain protein 7	
Immunogen Catalog Number: AG38516	Calculated MW: 55kDa, 489aa	
	Observed MW: 65-70 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : COLO 320 cells, RT-4 cells, U-251 cells
Species Specificity: human	

Background Information

UBXN7(UBX domain protein 7) is a protein with ubiquitin-binding activity and ubiquitin-proteasome-binding activity. It contains four different domains, namely, the N-terminal UBA domain, the UAS domain, the Ubiquitin Interaction motif (UIM) and the C-terminal UBX domain. The level of UBXN7 will affect the metabolic state of cells. High levels of UBXN7 and HIF-1 α accumulation support glycolysis, while inactivation of UBXN7 and activation of NRF2 increase oxidative phosphorylation (PMID: 33444648).

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 32722-1-AP (UBXN7 antibody) at dilution of 1:4000 incubated at room temperature for 1.5 hours.