

For Research Use Only

MBL Polyclonal antibody

Catalog Number: 32750-1-AP



Basic Information

Catalog Number:

32750-1-AP

Concentration:

400 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG3157

GenBank Accession Number:

NM_010776.1

GeneID (NCBI):

17195

UNIPROT ID:

P41317

Full Name:

mannose-binding lectin (protein C) 2

Calculated MW:

26 kDa

Observed MW:

26-30 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

mouse, rat

Positive Controls:

WB : mouse liver tissue, mouse serum, rat plasma

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 32750-1-AP (Mbl2 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.



Various lysates were subjected to SDS PAGE followed by western blot with 32750-1-AP (Mbl2 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.