

For Research Use Only

# RBM34 Polyclonal antibody

Catalog Number: 32834-1-AP



## Basic Information

Catalog Number:

32834-1-AP

Concentration:

700 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG37711

GenBank Accession Number:

BC029451

GeneID (NCBI):

23029

UNIPROT ID:

P42696

Full Name:

RNA binding motif protein 34

Calculated MW:

430 aa, 49 kDa

Observed MW:

55 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB: 1:5000-1:50000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HEK-293T cells, HeLa cells, Jurkat cells

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 32834-1-AP (RBM34 antibody) at dilution of 1:20000 incubated at room temperature for 1.5 hours.