

For Research Use Only

RBM34 Polyclonal antibody, PBS Only

Catalog Number: 32834-1-PBS



Basic Information

Catalog Number:

32834-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG37711

GenBank Accession Number:

BC029451

GeneID (NCBI):

23029

UNIPROT ID:

P42696

Full Name:

RNA binding motif protein 34

Calculated MW:

430 aa, 49 kDa

Observed MW:

55 kDa

Purification Method:

Antigen affinity Purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 32834-1-AP (RBM34 antibody) at dilution of 1:20000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 32834-1-PBS in a different storage buffer formulation.