

For Research Use Only

RPL37 Polyclonal antibody

Catalog Number:32865-1-AP



Basic Information

Catalog Number:

32865-1-AP

Concentration:

400 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34979

GenBank Accession Number:

NM_000997

GeneID (NCBI):

6167

UNIPROT ID:

P61927

Full Name:

ribosomal protein L37

Calculated MW:

11 kDa

Observed MW:

15-20 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB: 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HEK-293T cells, HeLa cells, Jurkat cells, Raji cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 32865-1-AP (RPL37 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.