

For Research Use Only

FAM173A Polyclonal antibody

Catalog Number: 32901-1-AP



Basic Information

Catalog Number:

32901-1-AP

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG37690

GenBank Accession Number:

BC001181

GeneID (NCBI):

65990

UNIPROT ID:

Q9BQD7

Full Name:

family with sequence similarity 173, member A

Calculated MW:

235 aa, 25 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

IHC: 1:50-1:500

Applications

Tested Applications:

IHC, ELISA

Species Specificity:

human, mouse

Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Positive Controls:

IHC: mouse kidney tissue,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded mouse kidney tissue slide using 32901-1-AP (FAM173A antibody) at dilution of 1:200 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunohistochemical analysis of paraffin-embedded mouse kidney tissue slide using 32901-1-AP (FAM173A antibody) at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).