

For Research Use Only

TBKBP1 Polyclonal antibody, PBS Only

Catalog Number: 32915-1-PBS



Basic Information

Catalog Number:

32915-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG38233

GenBank Accession Number:

NM_014726.2

GeneID (NCBI):

9755

UNIPROT ID:

A7MCY6

Full Name:

TBK1 binding protein 1

Calculated MW:

68kDa, 615aa

Observed MW:

80 kDa

Purification Method:

Antigen affinity Purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 32915-1-AP (TBKBP1 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 32915-1-PBS in a different storage buffer formulation.