

For Research Use Only

# EVI5L Polyclonal antibody, PBS Only

Catalog Number:32934-1-PBS



## Basic Information

Catalog Number:

32934-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG35996

GenBank Accession Number:

BC014111

GeneID (NCBI):

115704

UNIPROT ID:

Q96CN4

Full Name:

ecotropic viral integration site 5-like

Observed MW:

91,92 kDa

Purification Method:

Antigen affinity Purification

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

## Background Information

Evi5l (Ecotropic Viral Integration Site 5-Like), encodes a RabGAP that is known to block the formation of cilia when overexpressed (PMID: 26136363). EVI5L has ubiquitously low expression but is slightly higher in the testis and fetal brain.

## Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 32934-1-AP (EVI5L antibody) at dilution of 1:6000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 32934-1-PBS in a different storage buffer formulation.