

For Research Use Only

TTBK1 Polyclonal antibody, PBS Only

Catalog Number: 32938-1-PBS



Basic Information

Catalog Number: 32938-1-PBS	GenBank Accession Number: BC036764	Purification Method: Antigen affinity Purification
Concentration: 1 mg/ml	GeneID (NCBI): 84630	
Source: Rabbit	UNIPROT ID: Q5TCY1	
Isotype: IgG	Full Name: tau tubulin kinase 1	
Immunogen Catalog Number: AG36655	Calculated MW: 142/70 kDa	
	Observed MW: 200 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human, mouse, rat

Background Information

TTBK1 (Tau Tubulin kinase 1), also known as BDTK (brain-derived Tau kinase), is a 1,321 amino acid protein that contains one protein kinase domain and belongs to the serine/threonine protein kinase family. Localized to the cytoplasm and expressed at high levels in brain and at lower levels in testis and spinal cord, TTBK1 functions as a serine/threonine kinase that can phosphorylate Tau (a protein involved in Tubulin polymerization) on threonine, tyrosine and serine residues. Specifically, TTBK1 uses divalent cations, such as magnesium and manganese, to catalyze the ATP-dependent transfer of a phosphate group onto Tau, creating a phosphoprotein and ADP. Phosphorylation of Tau causes its aggregation and subsequent loss of function, suggesting an important role for TTBK1 in the control of Tubulin dynamics. Two isoforms of TTBK1 are expressed due to alternative splicing events. The molecular weight observed is consistent with what has been described in PMID:16923168.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 32938-1-AP (TTBK1 antibody) at dilution of 1:8000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 32938-1-PBS in a different storage buffer formulation.