

For Research Use Only

RNF185 Polyclonal antibody

Catalog Number:32956-1-AP



Basic Information

Catalog Number: 32956-1-AP	GenBank Accession Number: NM_152267.3	Purification Method: Antigen affinity Purification
Concentration: 400 ug/ml	GeneID (NCBI): 91445	Recommended Dilutions: WB: 1:500-1:2000
Source: Rabbit	UNIPROT ID: Q96GF1	
Isotype: IgG	Full Name: ring finger protein 185	
Immunogen Catalog Number: AG38937	Calculated MW: 20kDa, 192aa	
	Observed MW: 22 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : mouse testis tissue, rat testis tissue, human testis tissue
Species Specificity: human, mouse, rat	

Background Information

RNF185 is a RING finger domain-containing ubiquitin ligase implicated in ER-associated degradation. It is the MOM E3 ligase responsible for regulation of mitochondrial autophagy (PMID: 21931693).

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 32956-1-AP (RNF185 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.