

For Research Use Only

RNF185 Polyclonal antibody, PBS Only

Catalog Number:32956-1-PBS



Basic Information

Catalog Number: 32956-1-PBS	GenBank Accession Number: NM_152267.3	Purification Method: Antigen affinity Purification
Concentration: 1 mg/ml	GeneID (NCBI): 91445	
Source: Rabbit	UNIPROT ID: Q96GF1	
Isotype: IgG	Full Name: ring finger protein 185	
Immunogen Catalog Number: AG38937	Calculated MW: 20kDa, 192aa	
	Observed MW: 22 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human, mouse, rat

Background Information

RNF185 is a RING finger domain-containing ubiquitin ligase implicated in ER-associated degradation. It is the MOM E3 ligase responsible for regulation of mitochondrial autophagy (PMID: 21931693).

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 32956-1-AP (RNF185 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 32956-1-PBS in a different storage buffer formulation.