

For Research Use Only

FAM21C Polyclonal antibody

Catalog Number:33047-1-AP



Basic Information

Catalog Number:

33047-1-AP

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG15066

GenBank Accession Number:

BC006456

GeneID (NCBI):

253725

UNIPROT ID:

Q9Y4E1

Full Name:

family with sequence similarity 21, member C

Calculated MW:

1320 aa, 145 kDa

Observed MW:

200 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB: 1:1000-1:6000

IF/ICC: 1:200-1:800

Applications

Tested Applications:

WB, IF/ICC, ELISA

Species Specificity:

human

Positive Controls:

WB : HeLa cells, HepG2 cells, HuH-7 cells, L02 cells

IF/ICC : A431 cells,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 33047-1-AP (FAM21C antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.



Immunofluorescent analysis of (-20°C Methanol) fixed A431 cells using FAM21C antibody (33047-1-AP) at dilution of 1:400 and CoraLite® 488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2).