

For Research Use Only

NRG2 Polyclonal antibody

Catalog Number: 33064-1-AP



Basic Information

Catalog Number:

33064-1-AP

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG39953

GenBank Accession Number:

NM_004883.3

GeneID (NCBI):

9542

UNIPROT ID:

O14511-1

Full Name:

neuregulin 2

Calculated MW:

92 kDa

Observed MW:

46 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB: 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : RT-4 cells, U-251 cells, U-87 MG cells

Background Information

NRG2 is a novel member of the neuregulin family of growth and differentiation factors. Through interaction with the ERBB family of receptors, this protein induces the growth and differentiation of epithelial, neuronal, glial, and other types of cells. The gene consists of 12 exons and the genomic structure is similar to that of neuregulin 1, another member of the neuregulin family of ligands. The products of these genes mediate distinct biological processes by acting at different sites in tissues and eliciting different biological responses in cells. This gene is located close to the region for demyelinating Charcot-Marie-Tooth disease locus, but is not responsible for this disease.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 33064-1-AP (NRG2 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.