

For Research Use Only

Phosducin Polyclonal antibody

Catalog Number: 33132-1-AP



Basic Information

Catalog Number:

33132-1-AP

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG37137

GenBank Accession Number:

NM_002597.4

GeneID (NCBI):

5132

UNIPROT ID:

P20941

Full Name:

phosducin

Calculated MW:

28 kDa, 246 aa

Observed MW:

32 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB: 1:5000-1:50000

IHC: 1:50-1:500

Applications

Tested Applications:

WB, IHC, ELISA

Species Specificity:

human, mouse

Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Positive Controls:

WB : mouse retina tissue,

IHC : mouse eye tissue,

Background Information

Phosducin/PDC, which tightly binds betagamma-subunits of heterotrimeric G-proteins, has been conjectured to play a role in regulating second messenger signaling cascades. Phosducin may act as a phosphorylation-dependent switch in second messenger signaling cascades, regulating the kinetics of desensitization processes by controlling the activity of Gbetagamma-dependent GRKs (PMID: 9020189).

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 33132-1-AP (PDC antibody) at dilution of 1:20000 incubated at room temperature for 1.5 hours.



Immunohistochemical analysis of paraffin-embedded mouse eye tissue slide using 33132-1-AP (Phosducin antibody) at dilution of 1:200 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunohistochemical analysis of paraffin-embedded mouse eye tissue slide using 33132-1-AP (Phosducin antibody) at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).