

For Research Use Only

# EPS8L1 Polyclonal antibody, PBS Only

Catalog Number:33149-1-PBS



## Basic Information

Catalog Number:

33149-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG38144

GenBank Accession Number:

BC015763

GeneID (NCBI):

54869

UNIPROT ID:

Q8TE68

Full Name:

EPS8-Like 1

Calculated MW:

80 kDa

Observed MW:

80-85 kDa

Purification Method:

Antigen affinity Purification

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 33149-1-AP (EPS8L1 antibody) at dilution of 1:6000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 33149-1-PBS in a different storage buffer formulation.