

For Research Use Only

GTF3A Polyclonal antibody, PBS Only

Catalog Number: 33192-1-PBS



Basic Information

Catalog Number:

33192-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG35931

GenBank Accession Number:

NM_002097

GeneID (NCBI):

2971

UNIPROT ID:

Q92664

Full Name:

general transcription factor IIIA

Calculated MW:

41kDa

Observed MW:

42 kDa

Purification Method:

Antigen affinity Purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 33192-1-AP (GTF3A antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 33192-1-PBS in a different storage buffer formulation.