

For Research Use Only

# KCTD19 Polyclonal antibody

Catalog Number:33194-1-AP



## Basic Information

Catalog Number:

33194-1-AP

Concentration:

400 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG36660

GenBank Accession Number:

NM\_001100915.3

GeneID (NCBI):

146212

UNIPROT ID:

Q17RG1

Full Name:

potassium channel tetramerisation domain containing 19

Calculated MW:

105 kDa, 926aa

Observed MW:

105 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB: 1:1000-1:4000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : mouse testis tissue, rat testis tissue

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 33194-1-AP (KCTD19 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.