

For Research Use Only

MAFK Polyclonal antibody, PBS Only

Catalog Number: 33199-1-PBS



Basic Information

Catalog Number:

33199-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG38230

GenBank Accession Number:

NM_002360.3

GeneID (NCBI):

7975

UNIPROT ID:

O60675

Full Name:

v-maf musculoaponeurotic
fibrosarcoma oncogene homolog K
(avian)

Calculated MW:

18 kDa, 156aa

Observed MW:

18 kDa

Purification Method:

Antigen affinity Purification

Applications

Tested Applications:

WB, IF/ICC, Indirect ELISA

Species Specificity:

human

Background Information

The small Maf proteins (sMafs) are basic region leucine zipper (bZIP)-type transcription factors. MAFF, MAFK, and MAFG are the three vertebrate sMafs (PMID: 27058431). They act as transcriptional activators by dimerizing with other (usually larger) basic-zipper proteins, such as NFE2, NFE2L1/NRF1, NFE2L2/NRF2, and NFE2L3/NRF3, and recruiting them to specific DNA-binding sites (PMID: 8932385; 9150357). In cancer cell lines, MAFK promotes epithelial-to-mesenchymal transition (EMT) by downregulating E-cadherin expression through the targeting of transmembrane glycoprotein nonmetastatic B (GPNMB), a potent inducer of EMT (PMID: 28400538; 36061142).

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 33199-1-AP (MAFK antibody) at dilution of 1:4000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 33199-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed HeLa cells using MAFK antibody (33199-1-AP) at dilution of 1:800 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2), CL594-Phalloidin (red). This data was developed using the same antibody clone with 33199-1-PBS in a different storage buffer formulation.