

For Research Use Only

VSIG10 Polyclonal antibody, PBS Only

Catalog Number:33200-1-PBS



Basic Information

Catalog Number:

33200-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG37432

GenBank Accession Number:

BC034471

GeneID (NCBI):

54621

UNIPROT ID:

Q8NOZ9

Full Name:

hypothetical protein FLJ20674

Calculated MW:

60 kDa

Observed MW:

60 kDa

Purification Method:

Antigen affinity Purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, rat

Background Information

VSIG10 (V-set and immunoglobulin domain containing 10) is a protein-coding gene located on human chromosome 12q24.23 that encodes a transmembrane protein belonging to the immunoglobulin superfamily. VSIG10 is a protein whose expression is altered in response to SARS-CoV-2 infection, making it a novel candidate biomarker for COVID-19.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 33200-1-AP (VSIG10 antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 33200-1-PBS in a different storage buffer formulation.



Various lysates were subjected to SDS PAGE followed by western blot with 33200-1-AP (VSIG10 antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 33200-1-PBS in a different storage buffer formulation.