

For Research Use Only

# FAM21A Polyclonal antibody, PBS Only

Catalog Number:33208-1-PBS



## Basic Information

Catalog Number:

33208-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG36997

GenBank Accession Number:

BC082258

GeneID (NCBI):

387680

UNIPROT ID:

Q641Q2

Full Name:

family with sequence similarity 21,  
member A

Observed MW:

170 kDa

Purification Method:

Antigen affinity Purification

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 33208-1-AP (FAM21A antibody) at dilution of 1:8000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 33208-1-PBS in a different storage buffer formulation.