

For Research Use Only

DEFB119 Polyclonal antibody, PBS Only

Catalog Number: 33236-1-PBS



Basic Information

Catalog Number: 33236-1-PBS	GenBank Accession Number: BC062212	Purification Method: Antigen affinity Purification
Concentration: 1 mg/ml	GeneID (NCBI): 245932	
Source: Rabbit	UNIPROT ID: Q8N690	
Isotype: IgG	Full Name: defensin, beta 119	
Immunogen Catalog Number: AG36816	Calculated MW: 84 aa, 10 kDa	
	Observed MW: 10 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human

Background Information

DEFB119 is a reproductive tract-specific β -defensin highly expressed in the testis and the epididymis (16). Single-cell sequencing analysis of testicular cells isolated from nonobstructive azoospermia men demonstrates an elevated expression of DEFB119 in the Sertoli cells. Male mice lacking the mouse ortholog DEFB19 demonstrated subfertility in both sexes. DEFB119 is also expressed in sperm and is involved in sperm chemotaxis. Elevated seminal DEFB119 was associated with a decrease in the observed genera, diversity and evenness of bacterial communities, and further distortion of the metacommunities (PMID: 39359400).

Storage

Storage:
Store at -80°C .
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 33236-1-AP (DEFB119 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 33236-1-PBS in a different storage buffer formulation.