

For Research Use Only

ZNF70 Polyclonal antibody

Catalog Number:33264-1-AP



Basic Information

Catalog Number:

33264-1-AP

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG38011

GenBank Accession Number:

BC040161

GeneID (NCBI):

7621

UNIPROT ID:

Q9UC06

Full Name:

zinc finger protein 70

Observed MW:

50 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB: 1:1000-1:6000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB: Jurkat cells,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 33264-1-AP (ZNF70 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.