

For Research Use Only

# CDC2L2 Polyclonal antibody

Catalog Number:33277-1-AP



## Basic Information

Catalog Number:

33277-1-AP

Concentration:

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG36622

GenBank Accession Number:

NM\_001313896.1

GeneID (NCBI):

728642

UNIPROT ID:

Q9UQ88

Full Name:

cell division cycle 2-like 2 (PITSLRE proteins)

Calculated MW:

91kDa, 783aa

Purification Method:

Antigen affinity Purification

## Applications

Tested Applications:

ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data