

For Research Use Only

CDC2L2 Polyclonal antibody, PBS Only

Catalog Number:33277-1-PBS



Basic Information

Catalog Number:

33277-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG36622

GenBank Accession Number:

NM_001313896.1

GeneID (NCBI):

728642

UNIPROT ID:

Q9UQ88

Full Name:

cell division cycle 2-like 2 (PITSLRE proteins)

Calculated MW:

91kDa, 783aa

Purification Method:

Antigen affinity Purification

Applications

Tested Applications:

Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data