

For Research Use Only

MAG/Siglec-4a Polyclonal antibody

Catalog Number:33295-1-AP



Basic Information

Catalog Number:

33295-1-AP

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG3251

GenBank Accession Number:

NM_017190.4

GeneID (NCBI):

29409

UNIPROT ID:

P07722-1

Full Name:

myelin-associated glycoprotein

Calculated MW:

69 kDa

Observed MW:

70-100 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB: 1:2000-1:10000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

mouse, rat

Positive Controls:

WB : mouse brain tissue, rat brain tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 33295-1-AP (MAG/Siglec-4a antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.