

For Research Use Only

C9orf169 Polyclonal antibody, PBS Only

Catalog Number: 33316-1-PBS



Basic Information

Catalog Number: 33316-1-PBS	GenBank Accession Number: BC052297	Purification Method: Antigen affinity Purification
Concentration: 1 mg/ml	GeneID (NCBI): 375791	
Source: Rabbit	UNIPROT ID: A8MQ03	
Isotype: IgG	Full Name: chromosome 9 open reading frame 169	
Immunogen Catalog Number: AG37408	Calculated MW: 144 aa, 15 kDa	
	Observed MW: 18 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human

Background Information

CYSRT1 is highly conserved among vertebrates and has putative antimicrobial activity. CYSRT1 colocalizes with LCE proteins and contributes to the constitutive epidermal antimicrobial host defense repertoire (PMID: 36804407).

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 33316-1-AP (C9orf169 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 33316-1-PBS in a different storage buffer formulation.