

For Research Use Only

TMCO2 Polyclonal antibody

Catalog Number:33338-1-AP



Basic Information

Catalog Number:

33338-1-AP

Concentration:

400 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG37348

GenBank Accession Number:

NM_001008740.3

GeneID (NCBI):

127391

UNIPROT ID:

Q7Z6W1

Full Name:

transmembrane and coiled-coil domains 2

Calculated MW:

20kDa, 182aa

Observed MW:

20-25 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB: 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : mouse testis tissue, rat testis tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 33338-1-AP (TMCO2 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.