

For Research Use Only

TMCO2 Polyclonal antibody, PBS Only

Catalog Number:33338-1-PBS



Basic Information

Catalog Number:

33338-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG37348

GenBank Accession Number:

NM_001008740.3

GeneID (NCBI):

127391

UNIPROT ID:

Q7Z6W1

Full Name:

transmembrane and coiled-coil domains 2

Calculated MW:

20kDa, 182aa

Observed MW:

20-25 kDa

Purification Method:

Antigen affinity Purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 33338-1-AP (TMC02 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 33338-1-PBS in a different storage buffer formulation.