

For Research Use Only

# FCRL2 Polyclonal antibody

Catalog Number: 33416-1-AP



## Basic Information

Catalog Number:

33416-1-AP

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG4773

GenBank Accession Number:

NM\_030764.3

GeneID (NCBI):

79368

UNIPROT ID:

Q96LA5-1

Full Name:

Fc receptor-like 2

Calculated MW:

56 kDa

Observed MW:

60-70 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB: 1:500-1:2000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HeLa cells,

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 33416-1-AP (FCRL2 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.