

For Research Use Only

VSTM2B Polyclonal antibody

Catalog Number: 33423-1-AP



Basic Information

Catalog Number: 33423-1-AP	GenBank Accession Number: NM_001146339.1	Purification Method: Antigen affinity Purification
Source: Rabbit	GeneID (NCBI): 342865	Recommended Dilutions: WB: 1:500-1:2000
Isotype: IgG	UNIPROT ID: A6NLU5	
Immunogen Catalog Number: EG6734	Full Name: V-set and transmembrane domain containing 2B	
	Calculated MW: 30 kDa	
	Observed MW: 30 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : fetal human brain tissue,
Species Specificity: human	

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Fetal human brain tissue were subjected to SDS PAGE followed by western blot with 33423-1-AP (VSTM2B antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.