

For Research Use Only

ZNF85 Polyclonal antibody

Catalog Number: 33448-1-AP



Basic Information

Catalog Number:

33448-1-AP

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG38326

GenBank Accession Number:

NM_003429.4

GeneID (NCBI):

7639

UNIPROT ID:

Q03923

Full Name:

zinc finger protein 85

Calculated MW:

69kDa, 595aa

Observed MW:

70 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB: 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : HuH-7 cells, Daudi cells, Huh7 cells, SKOV-3 cells, mouse brain tissue, rat brain tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 33448-1-AP (ZNF85 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.



Various lysates were subjected to SDS PAGE followed by western blot with 33448-1-AP (ZNF85 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.