

For Research Use Only

SCAI Polyclonal antibody

Catalog Number: 33449-1-AP



Basic Information

Catalog Number:

33449-1-AP

Concentration:

300 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG38169

GenBank Accession Number:

NM_001144877.2

GeneID (NCBI):

286205

UNIPROT ID:

Q8N9R8

Full Name:

chromosome 9 open reading frame

126

Calculated MW:

70kDa, 606aa

Observed MW:

60 kDa

Purification Method:

Antigen affinity Purification

Recommended Dilutions:

WB: 1:1000-1:8000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : Jurkat cells, mouse brain tissue, rat brain tissue

Background Information

SCAI (suppressor of cancer cell invasion) is a conserved protein that regulates invasive cell migration through three-dimensional matrices. SCAI acts on the RhoA-Dia1 signal transduction pathway and localizes in the nucleus, where it binds and inhibits the myocardin-related transcription factor MAL by forming a ternary complex with serum response factor (SRF). Genome-wide expression analysis surprisingly reveals that one of the strongest upregulated genes after suppression of SCAI is beta1-integrin. Decreased levels of SCAI are tightly correlated with increased invasive cell migration, and SCAI is downregulated in several human tumours. Functional analysis of the beta1-integrin gene strongly argues that SCAI is a novel transcriptional cofactor that controls gene expression downstream of Dia1 to dictate changes in cell invasive behaviour (PMID: 19350017).

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 33449-1-AP (SCA1 antibody) at dilution of 1:4000 incubated at room temperature for 1.5 hours.