

For Research Use Only

# CCDC58 Polyclonal antibody, PBS Only

Catalog Number:33455-1-PBS



## Basic Information

Catalog Number:

33455-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG38762

GenBank Accession Number:

BC062725

GeneID (NCBI):

131076

UNIPROT ID:

Q4VC31

Full Name:

coiled-coil domain containing 58

Observed MW:

16 kDa

Purification Method:

Antigen affinity Purification

## Applications

Tested Applications:

WB, IP, Indirect ELISA

Species Specificity:

human, mouse, rat

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 33455-1-AP (CCDC58 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 33455-1-PBS in a different storage buffer formulation.



IP result of anti-CCDC58 (IP:33455-1-AP, 4ug; Detection:33455-1-AP 1:1000) with rat kidney tissue lysate 3160 ug. This data was developed using the same antibody clone with 33455-1-PBS in a different storage buffer formulation.